Aesculap Orthopaedics

Plasmapore® μ -CaP

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Inhalt

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- (1) **Volkman R, Eingartner C, Winter E, Weise K, Weller S**
Midterm results in 500 titanium alloy straight femoral stem prostheses – cemented and cementless technique
Eur J Orthop Surg Traumatol 8 (1998) 133–139
- (2) **Asmuth T, Bachmann J, Eingartner C, Feldmann C, Fünten aus der K, Holz F, Hübenthal L, Papp J, Quack G, Sauer G**
Results with the cementless Bicontact stem – multicenter study of 553 cases in: S. Weller, R. Volkman (editors) *The Bicontact hip system*
Georg Thieme Stuttgart (1998) 63–74
- (3) **Eingartner C, Volkman R, Winter E, Maurer F, Sauer G, Weller S, Weise K**
Results of an uncemented straight femoral shaft prosthesis after 9 years follow-up
J Arthroplasty 4 (2000) 440–447
- (4) **Badhe NP, Quinell RC, Howard PW**
The Uncemented Bi-Contact Total Hip Arthroplasty
J Arthroplasty 7 (2002) 896–901
- (5) **Gordon A, Wilkinson JM, Stockley I**
Midterm Outcome of the Plasmacup in Total Hip Arthroplasty
Hip International 12:2 (2002) 119–125
- (6) **Winkler-Gniewek W, Stallforth H, Ungethüm M**
Die Plasmapore-Beschichtung von Gelenkendothesen – ein neues Konzept. in: G. Friedebold (Editor) *Oberflächenstrukturierte Prothesen aus technischer und medizinischer Sicht*.
DVM, Berlin (1987) 59–70
- (7) **Gross U, Müller-Mai C, Fritz T, Voigt C, Knarse W, Schmitz HJ**
Implant surface roughness and mode of load transmission influence periimplant bone structure
in: G. Heimke, U. Soltesz and AJC Lee, Elsevier (editors) *Clinical Implant Materials, Advances in Biomaterials* 9, Amsterdam (1990) 303–308
- (8) **Fink U**
Plasmapore: a plasma-sprayed microporous titanium coating to improve the long term stability
in: D. Mainard (editor) *Actualités en Biomateriaux Vol. III*
Romillat Edit., Paris (1996) 97–104
- (9) **Bobyn JD, Pilliar RM, Cameron HU, Weatherly GC**
The optimum pore size for the fixation of porous surfaced metal implants by the ingrowth of bone
Clin Orthop 150 (1980) 263–270
- (10) **Pilliar RM, Lee JM, Maniopoulos C**
Observations on the effect of movement on bone ingrowth into porous-surfaced implants
Clin Orthop 208 (1986) 108–113
- (11) **Turner TM, Sumner DR, Urban RM, Rivero DP, Galante JO**
A comparative study of porous coatings in a weight-bearing total hip-arthroplasty model
J Bone Joint Surg A 68:9 (1986) 1396–1409
- (12) **Pitto RP, Böhner J, Hofmeister V**
Einflußgrößen der Primärstabilität acetabulärer Komponenten. Eine In-vitro-Studie
Biomed Technik 42 (1997) 363 – 368
- (13) **Weller S, Braun A, Gellrich JC, Gross U**
Importance of prosthesis design and surface structure for primary and secondary stability of uncemented hip joint prostheses
in: I.D. Learmonth (editor) *Interfaces in Total Hip Arthroplasty*
Springer London (1999) 81–101
- (14) **Jaffe WL, Scott DF**
Total hip arthroplasty with hydroxyapatite-coated prostheses
J Bone Joint Surg A 78:12 (1996) 1918–1934
- (15) **Petit R**
The use of hydroxyapatite in orthopaedic surgery: a ten-year review
Eur J Orthop Surg Traumatol 9 (1999) 71–74
- (16) **Hardy DCR, Frayssinet P, Delince PE**
Osteointegration of hydroxyapatite-coated stems of femoral prostheses
Eur J Orthop Surg Traumatol 9 (1999) 75–81
- (17) **McNally SA, Shepperd JA, Mann CV, Walczak JP**
The results at nine to twelve years of the use of a hydroxyapatite-coated femoral stem
J Bone Joint Surg B 82:3 (2000) 378–382
- (18) **Furlong RJ, Osborn JF**
Fixation of hip prostheses by hydroxyapatite ceramic coatings
J Bone Joint Surg B 73:5 (1991) 741–745
- (19) **Geesink RG, Hoefnagels NH**
Six-year results of hydroxyapatite-coated total hip replacement
J Bone Joint Surg B 77:4 (1995) 534–547
- (20) **Capello WN, D'Antonio JA, Feinberg JR, Manley MT**
Hydroxyapatite-coated femoral components in patients less than fifty years old. Clinical and radiographic results after five to eight years of follow-up
J Bone Joint Surg A 79:9 (1997) 1023–1029
- (21) **Vidalain JP, and the ARTRO group**
The Corail system in primary THA: results, lessons and comments from the series performed by the ARTRO group
Eur J Orthop Surg Traumatol 9 (1999) 87–90
- (22) **Epinette JA**
HA-coated hip implants: a 10 year follow-up
Eur J Orthop Surg Traumatol 9 (1999) 83–85
- (23) **Cook SD, Thomas KA, Kay SF, Jarcho M**
Hydroxyapatite-coated porous titanium for use as an orthopaedic biologic attachment system
Clin Orthop 230 (1988) 303–312
- (24) **Willmann G**
Coating of implants with hydroxyapatite-material connections between bone and metal
Adv Eng Mat 1:2 (1999) 95–105
- (25) **Karabatsos B, Myerthall SL, Fornasier VL, Binnington A, Maistrelli GL**
Osseointegration of hydroxyapatite porous-coated femoral implants in a canine model
Clin Orthop 392 (2001) 442–449
- (26) **Ishizawa H, Fujino M, Ogino M**
Histomorphometric evaluation of the thin hydroxyapatite layer formed through anodization followed by hydrothermal treatment
J Biomed Mater Res. 35:2 (1997) 199–206
- (27) **Wheeler DL, Capbell AA, Graff GL, Miller GJ**
Histological and biomechanical evaluation of calcium phosphate coatings applied through surface-induced mineralization to porous titanium implants
J Biomed Mater Res 34:4 (1997) 539–543
- (28) **Overgaard S, Bromose U, Lind M, Bunker C, Soballe K**
The influence of crystallinity of the hydroxyapatite coating on the fixation of implants. Mechanical and histomorphometric results
J Bone Joint Surg B 81:4 (1999) 725–731
- (29) **Thanner J, Kärrholm J, Herberts P, Malchau H**
Hydroxyapatite and tricalcium phosphate-coated cups with and without screw fixation
J Arthroplasty 15:4 (2000) 405–412
- (30) **Ducheyne P, Beight J, Cuckler J, Evans B, Radin S**
Effect of calcium phosphate coating characteristics on early post-operative bone tissue ingrowth
Biomaterials 11:8 (1990) 531–40

- (31) Bloebaum RD, Bachus KN, Rubman MH, Dorr LD
Postmortem comparative analysis of titanium and hydroxyapatite porous-coated femoral implants retrieved from the same patient
J Arthroplasty 8:2 (1993) 203-211
- (32) Moroni A, Caja VL, Egger EL, Trinchese L, Chao EYS
Histomorphometry of hydroxyapatite coated and uncoated porous titanium bone implants
Biomaterials 15:11 (1994) 926-930
- (33) Szmukler-Moncler S, Reingewirtz Y, Weber HP
Bone response to early loading: The effect of surface state,
in: Z. Davidovitch and J. Mah (editors) Biological mechanisms of tooth eruption, resorption and replacement by implants, Harvard society for the advancements of orthodontics, Boston (1998) 611-616
- (34) Collier JP, Surprenant VA, Mayor MB, Wrona M, Jensen RE, Surprenant HE
Loss of hydroxyapatite coating on retrieved total hip components
J Arthroplasty 8:4 (1993) 389-393
- (35) Bloebaum RD, Dupont JA
Osteolysis from a press-fit hydroxyapatite-coated implant
J Arthroplasty 8:2 (1993) 195-202
- (36) Bloebaum RD, Beeks D, Dorr LD, Savory CG, DuPont JA, Hofmann AA
Complications with hydroxyapatite particulate separation in total hip arthroplasty
Clin Orthop 298 (1994) 19-26
- (37) David A, Eitenmüller J, Muhr G, Pommer A, Bör HF, Ostermann PAW, Schildhauer TA
Mechanical and histological evaluation of hydroxyapatite-coated, titanium-coated and grit-blasted surfaces under weight bearing conditions
Arch Orthop Trauma Surg 114 (1995) 112-118
- (38) Furlong R
Severe osteolysis after third-body wear due to hydroxyapatite particles from acetabular cup coating
J Bone Joint Surg B 80:4 (1998) 745-746
- (39) Manley MT, Capello WN, D'Antonio JA, Edidin AA, Geesink RG
Fixation of acetabular cups without cement in total hip arthroplasty. A comparison of three different implant surfaces at a minimum duration of follow-up of five years
J Bone Joint Surg A 80:8 (1998) 1175-85
- (40) Galois L, Mainard D, Bordji K, Membre H, Marchal L, Foliguet B, Clement D, Delagoutte JP
Influence de la taille des pores sur la réhabilitation osseuse de 2 céramiques phospho-calciques: l'hydroxyapatite et le phosphate tricalcique
in: D. Mainard (editor) Actualités en Biomatériaux Vol. III
Romillat Edit., Paris (1996) 361-380
- (41) Cheang P, Khor KA
Addressing processing problems associated with plasma spraying of hydroxyapatite coatings
Biomaterials 17:5 (1996) 537-544

- (42) Park E, Condrate RA, Hoelzer DT
Interfacial characterization of plasma-spray coated calcium phosphate on Ti-6Al-4V
Journal of material Science: Materials in medicine 9 (1998) 643-649
- (43) Barrere F, Van Der Valk CM, Dalmeijer RA, Meijer G, Van Blitterswijk CA, De Groot K, Layrolle P
Osteogenicity of octacalcium phosphate coatings applied on porous metal implants
J Biomed Mater Res. 66A:4 (2003) 779-788
- (44) Chow LC, Markovic M, Takagi S
A dual constant-composition titration system as an in vitro resorption model for comparing dissolution rates of calcium phosphate biomaterials
J Biomed Mater Res. 65B:2 (2003) 245-251
- (45) Kumar M, Dasarathy H, Riley C
Electrodeposition of brushite coatings and their transformation to hydroxyapatite in aqueous solutions
J Biomed Mater Res 45:4 (1999) 302-310
- (46) Szmukler-Moncler S, Perrin D, Piattelli A, Scarano A
Evaluation of a soluble calcium phosphate coating obtained by electrochemical deposition: a pilot study in the pig maxillae
in: Z. Davidovitch and J. Mah (editors) Biological mechanisms of tooth eruption, resorption and replacement by implants, Harvard society for the advancements of orthodontics, Boston (1998) 481-485
- (47) Szmukler-Moncler S, Perrin D, Ahossi V, Pointaire Ph
Evaluation of a resorbable Calcium Phosphate (CaP) coating obtained by electrochemical deposition after 6 weeks of healing. A pilot study in the pig maxilla
Poster presentation 15th European Conference for Biomaterials, Bordeaux September 1999
- (48) Redepenning J, Schlessinger T, Burnham S, Lippiello L, Miyano J
Characterization of electrolytically prepared brushite and hydroxyapatite coatings on orthopedic alloys
J Biomed Mater Res 30:3 (1996) 287-294
- (49) Caulier H, van der Waerden JP, Wolke JG, Kalk W, Naert I, Jansen JA
A histological and histomorphometrical evaluation of the application of screw-designed calciumphosphate (Ca-P)-coated implants in the cancellous maxillary bone of the goat
J Biomed Mater Res 35:1 (1997) 19-30
- (50) Gottlander M, Johansson CB, Wennerberg A, Albrektsson T, Radin S, Ducheyne P
Bone tissue reactions to an electrophoretically applied calcium phosphate coating
Biomaterials 18:7 (1997) 551-557

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Aesculap AG & Co. KG