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**AESCULAP® PROSPACE® 3D**  
ZWISCHENWIRBELKÖRPERIMPLANTAT ZUR  
POSTERIOREN THORAKOLUMBALEN FUSION  
OPERATIONSTECHNIK

# AESCULAP® PROSPACE® 3D

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# AESCULAP® PROSPACE® 3D

## REFERENZEN

- (1) Tropiano P, Bronsard JJ, Louis C, Tallet JM, Sauget Y. Three column stabilisation through a posterior approach with a titanium PLASMAPORE® intervertebral block (PROSPACE®). Radiological and clinical study after 4 years. *Rivista di Neuroradiologia*. 1999;12(Suppl 1):89-94.
- (2) Kroppenstedt S, Gulde M, Schönmayr R. Radiological comparison of instrumented posterior lumbar interbody fusion with one or two closed-box PLASMAPORE® coated titanium cages. Follow-up study over more than seven years. *Spine*. 2008;33(19):2083-8.
- (3) Kreinest M, Schmahl D, Grützner PA, Matschke S. Radiological Results and Clinical Patient Outcome After Implantation of a Hydraulic Expandable Vertebral Body Replacement following Traumatic Vertebral Fractures in the Thoracic and Lumbar Spine: A 3-Year Follow-Up. *Spine (Phila Pa 1976)*. 2017 Apr 15;42(8):E482-E489.
- (4) Takeuchi M, Yasuda M, Niwa A, Wakao N, Nakura T, Osuka K et al. PLASMAPORE®-coated titanium cervical cages induce more rapid and complete bone fusion after anterior cervical discectomy and fusion as compared to noncoated titanium cage. *World Neurosurgery*. 2014;82(3/4):519-22.
- (5) Vanek P, Bradac O, Konopkova R, de Lacy P, Lacman J, Benes V. Treatment of thoracolumbar trauma by short-segment percutaneous transpedicular screw instrumentation: prospective comparative study with a minimum 2-year follow-up. *J Neurosurg Spine*. 2014;20:150-6.
- (6) Beisse R. Endoscopic surgery on the thoracolumbar junction of the spine. *Eur Spine J*. 2006;15:687-704.
- (7) Finger T, Bayerl S, Onken J, Czabanka M, Woitzik J, Vajkoczy P. Sacropelvic fixation versus fusion to the sacrum for spondylodesis in multilevel degenerative spine disease. *Eur Spine J*. 2014;23:1013-20.
- (8) Ngoc Bich VU et al. In vitro and in vivo biocompatibility of Ti-6Al-4V titanium alloy and UHMWPE polymer for total hip replacement. *Biomedical Research and Therapy*. 2016;3(3):567-77.
- (9) Van der Stok J, Van der Jagt O, Yavari S, De Haas M, Waarsing J, Jahr H et al. Selective laser melting-produced porous titanium scaffolds regenerate bone in critical size cortical bone defects. *Journal of Orthopaedic Research*. 2013;31(5):792-9.
- (10) In vivo study of porous metallic lattice structures, Ulm, 2011. The biocompatibility, osseointegration and biomechanical properties of porous Ti6Al4V implants manufactured by SLM were tested under mechanical loading conditions in an ovine model study sponsored by Aesculap AG. The samples were evaluated histologically and biomechanically after implantation. Porous Ti6Al4V implants exhibited very good biocompatibility, bone-implant interface strength and osseointegration. Six months after implantation, bone ingrowth on and into the porous Ti6Al4V implants was reported. Inflammatory reactions that may influence bone formation were not observed.
- (11) Karageorgiou V, Kaplan D. Porosity of 3D biomaterial scaffolds and osteogenesis. *Biomaterials*. 2005;26(27):5474-91.
- (12) Taniguchi N et al. Effect of pore size on bone ingrowth into porous titanium implants fabricated by additive manufacturing: an in vivo experiment. *Materials Science and Engineering: C*. 2016;59:690-701.
- (13) Van Bael S et al. The effect of pore geometry on the in vitro biological behavior of human periosteum-derived cells seeded on selective laser-melted Ti6Al4V bone scaffolds. *Acta biomaterialia*. 2012;8(7):2824-34.
- (14) Elias CN et al. Mechanical and clinical properties of titanium and titanium-based alloys (Ti G2, Ti G4 cold worked nano structured and Ti G5) for biomedical applications. *Journal of Materials Research and Technology*. 2019;8(1):1060-9.
- (15) Kuhn JL, Goldstein SA, Choi K, London M, Feldkamp LA, Matthews LS. Comparison of the trabecular and cortical tissue moduli from human iliac crests. *J Orthop Res*. 1989;7(6):876-84.

- (16) Ratner BD, Hoffmann AS, Schoen FJ, Lemons JE. An Introduction to Materials in Medicine. Academic Press. 1996.
- (17) Internal Report U-0464/19: Product Qualification Validation of Aesculap 3D Cages lattice structure, Tuttlingen, 2019. In the scope of a product-related PQ validation, test specimens were manufactured with a lattice structure corresponding to the Aesculap 3D Cages. Subsequently, the strength was determined by a static compression test, the microstructure by a structural analysis as well as the dimensional accuracy by measuring the lattice struts. The investigated lattice test specimens totally satisfied the acceptance criteria in terms of lattice dimensions, strength and microstructure.
- (18) Chen Y, Wang X, Lu X, Yang L, Yang H, Yuan W et al. Comparison of titanium and polyetheretherketone (PEEK) cages in the surgical treatment of multilevel cervical spondylotic myelopathy: a prospective, randomized, control study with over 7-year follow-up. *Eur Spine J.* 2013;22(7):1539–46.
- (19) Brizuela A et al. Influence of the elastic modulus on the Osseointegration of Dental Implants. *Materials.* 2019;12(6):980.
- (20) Azami M, Moztarzadeh F, Tahriri M. Preparation, characterization and mechanical properties of controlled porous gelatin/hydroxyapatite nanocomposite through layer solvent casting combined with freeze-drying and lamination techniques. *Journal of Porous Materials.* 2010;17(3):313–20.
- (21) Rae PJ, Brown EN, Orlor EB. The mechanical properties of poly (ether-ether-ketone) (PEEK) with emphasis on the large compressive strain response. *Polymer.* 2007;48(2):598–615.
- (22) Internal Report U-0246/19, Tuttlingen 2019: Equivalence test on 3D Cages, compressive specimens with 3D Cage lattice structure and standard tensile specimens. The aim of the study was to verify the equivalence of mechanical tests on both test specimens and Aesculap 3D Cages. For this purpose, test specimens and implants were built up in a single production job by using laser sintering technology and were tested afterwards. The results of the static and dynamic compression tests on lattice test specimens as well as tensile tests on tensile specimens showed a high equivalence to the results of the "Dynamic Axial Compression Test according to ASTM F2077" on implants. The test specimens showed a well-known behavior, and all nominal or standard values were reached or exceeded.
- (23) Usability-Test, Usability Validation of AESCULAP® PROSPACE® 3D Cages, Tübingen, 2019. The usability of the AESCULAP® 3D Cage System PROSPACE® 3D was tested in April 2019, in a cadaver workshop with six independent test persons as intended users (surgeons specialized in spinal surgery or comparable fields). Parameters such as implant visibility under x-ray control, mechanical stability of the implant/instrument interface and implant surface evaluation in terms of tissue injury risk were tested among others. Acceptance criteria were fulfilled for all the above-mentioned parameters. All test users confirmed the absence of critical features that must be improved prior to clinical use.
- (24) Rehnitz, Christoph, PD Dr. med. Radiological image evaluation of AESCULAP® interbody fusion devices, Heidelberg, 2019. CT and X-ray visualization of different AESCULAP® interbody fusion cages (full titanium, porous Ti6Al4V and PLASMAPORE® XP cages) was tested in a cadaver setup. A radiologist evaluated the implant visibility and the presence of artefacts that may limit the visualization of adjacent structures. Visualization and assessment of implant position was achieved in X-ray and CT for all tested cages. Minor artefacts were visible in CT reconstructions in the surrounding of porous Ti6Al4V and full titanium implants. Porous Ti6Al4V implants showed slightly fewer artefacts in CT in comparison to full titanium implants. The minor artefacts observed did not limit the assessment of the surrounding anatomical structures.

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[www.facebook.com/bbraunat](https://www.facebook.com/bbraunat) | [www.instagram.com/bbraun\\_at](https://www.instagram.com/bbraun_at) | [www.linkedin.com/company/bbraun-group](https://www.linkedin.com/company/bbraun-group)

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B. Braun Austria GmbH | Otto Braun-Straße 3-5 | 2344 Maria Enzersdorf
Hersteller nach MDR ist Aesculap AG, Am Aesculap-Platz, 78532 Tuttlingen, Deutschland