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# AESCULAP® CeSPACE® 3D

## REFERENZEN

- (1) Suchomel P, Jurák L, Antinheimo J, Pohjola J, Stulik J, Meisel HJ, et al. Does sagittal position of the CTDR-related centre of rotation influence functional outcome? Prospective 2-year follow-up analysis. *Eur Spine J.* 2014 May;23(5):1124-34.
- (2) Meisel HJ, Jurák L, Antinheimo J, Arregui R, Bruchmann B, Čabreja M, et al. Four-year results of a prospective single-arm study on 200 semi-constrained total cervical disc prostheses: clinical and radiographic outcome. *J Neurosurg Spine.* 2016 Jun;3:1-10.
- (3) Boselie TFM, van Santbrink H, de Bie RA, van Mameren H. Pilot Study of Sequence of Segmental Contributions in the Lower Cervical Spine During Active Extension and Flexion: Healthy Controls Versus Cervical Degenerative Disc Disease Patients. *Spine (Phila Pa 1976).* 2017 Jun 1;42(11):E642-E647.
- (4) Epstein NE. Anterior cervical dynamic ABC plating with single level corpectomy and fusion in forty-two patients. *Spinal Cord.* 2003 Mar;41(3):153-8. PubMed PMID:12612617.
- (5) Stulik J, Pitzen TR, Chrobok J, Ruffing S, Drumm J, Sova L, et al. Fusion and failure following anterior cervical plating with dynamic or rigid plates: 6-months results of a multicentric, prospective, randomized, controlled study. *Euro Spine J.* 2007;16:1689-94.
- (6) Krayenbühl N, Schneider C, Landolt H, Fandino J. Use of an empty, Plasmapore-covered titanium cage for interbody fusion after anterior cervical microdiscectomy. *J Clin Neurosci.* 2008; 15(1):11-7.
- (7) Takeuchi M, Yasuda M, Niwa A, Wakao N, Nakura T, Osuka K, et al. Plasmapore-coated titanium cervical cages induce more rapid and complete bone fusion after anterior cervical discectomy and fusion as compared to noncoated titanium cage. *World Neurosurg.* 2014;82(3-4):519-22.
- (8) Stulik J, Vyskocil T, Sebesta P, Kryl J. Atlantoaxial fixation using the polyaxial screw-rod system. *Eur Spine J.* 2007 Apr; 16(4):479-84. PubMed PMID: 17051397; PubMed Central PMCID: PMC2229812.
- (9) Tippets RH, Apfelbaum RI. Anterior cervical fusion with the Caspar instrumentation system. *Neurosurgery.* 1988 Jun; 22(6 Pt 1):1008-13. PubMed PMID: 3419561.
- (10) Stulik J, Nesnidal P, Kryl J, Vyskocil T, Barna M. Kyphotic deformities of the cervical spine. 28th Annual Meeting of the AANS/CNS Section on Disorders of the cervical Spine and peripheral Nerves. March 2012 Orlando, Florida.
- (11) Ngoc Bich VU, et al. In vitro and in vivo biocompatibility of Ti-6Al-4V titanium alloy and UHMWPE polymer for total hip replacement. *Biomedical Research and Therapy.* 2016;3(3):567-77.
- (12) Van der Stok J, Van der Jagt O, Yavari S, De Haas M, Waarsing J, Jahr H, et al. Selective laser melting-produced porous titanium scaffolds regenerate bone in critical size cortical bone defects. *Journal of Orthopaedic Research* 2013;31(5):792-9.
- (13) In vivo study of porous metallic lattice structures, Ulm, 2011. The biocompatibility, osseointegration and biomechanical properties of porous Ti6Al4V implants manufactured by SLM were tested under mechanical loading conditions in an ovine model study sponsored by Aesculap AG. The samples were evaluated histologically and biomechanically after implantation. Porous Ti6Al4V implants exhibited very good biocompatibility, bone-implant interface strength and osseointegration. Six months after implantation, bone ingrowth on and into the porous Ti6Al4V implants was reported. Inflammatory reactions that may influence bone formation were not observed.
- (14) Karageorgiou V, Kaplan D. Porosity of 3D biomaterial scaffolds and osteogenesis. *Biomaterials.* 2005;26(27):5474-91.
- (15) Taniguchi N, et al. Effect of pore size on bone ingrowth into porous titanium implants fabricated by additive manufacturing: an in vivo experiment. *Materials Science and Engineering: C.* 2016;59:690-701.
- (16) van Bael S, et al. The effect of pore geometry on the in vitro biological behavior of human periosteum-derived cells seeded on selective laser-melted Ti6Al4V bone scaffolds. *Acta biomaterialia.* 2012;8(7):2824-34.

- (17) Elias CN, et al. Mechanical and clinical properties of titanium and titanium-based alloys (Ti G2, Ti G4 cold worked nano-structured and Ti G5) for biomedical applications. *Journal of Materials Research and Technology*. 2019;8(1):1060-9.
- (18) Kuhn JL, Goldstein SA, Choi K, London M, Feldkamp LA, Matthews LS. Comparison of the trabecular and cortical tissue moduli from human iliac crests. *J Orthop Res*. 1989;7(6):876-84.
- (19) Ratner BD, Hoffmann AS, Schoen FJ, Lemons JE. *An Introduction to Materials in Medicine*. Academic Press. 1996
- (20) Internal Report U-0464/19: Product Qualification Validation of Aesculap 3D Cages lattice structure, Tuttlingen, 2019. In the scope of a product-related PQ validation, test specimens were manufactured with a lattice structure corresponding to the Aesculap 3D Cages. Subsequently, the strength was determined by a static compression test, the microstructure by a structural analysis as well as the dimensional accuracy by measuring the lattice struts. The investigated lattice test specimens totally satisfied the acceptance criteria in terms of lattice dimensions, strength and microstructure.
- (21) Chen Y, Wang X, Lu X, Yang L, Yang H, Yuan W, et al. Comparison of titanium and polyetheretherketone (PEEK) cages in the surgical treatment of multilevel cervical spondylotic myelopathy: a prospective, randomized, control study with over 7-year follow-up. *Eur Spine J*. 2013;22(7):1539-46.
- (22) Brizuela A, et al. Influence of the elastic modulus on the Osseointegration of Dental Implants. *Materials*. 2019;12(6):980.
- (23) Azami M, Moztafzadeh F, Tahriri M. Preparation, characterization and mechanical properties of controlled porous gelatin/hydroxyapatite nanocomposite through layer solvent casting combined with freeze-drying and lamination techniques. *Journal of Porous Materials*. 2010;17(3):313-20.
- (24) Rae PJ, Brown EN, Orler EB. The mechanical properties of poly (ether-ether-ketone) (PEEK) with emphasis on the large compressive strain response. *Polymer*. 2007;48(2):598-615.
- (25) Internal Report U-0246/19, Tuttlingen 2019: Equivalence test on 3D Cages, compressive specimens with 3D Cage lattice structure and standard tensile specimens. The aim of the study was to verify the equivalence of mechanical tests on both test specimens and Aesculap 3D Cages. For this purpose, test specimens and implants were built up in a single production job by using laser sintering technology and were tested afterwards. The results of the static and dynamic compression tests on lattice test specimens as well as tensile tests on tensile specimens showed a high equivalence to the results of the "Dynamic Axial Compression Test according to ASTM F2077" on implants. The test specimens showed a well-known behavior, and all nominal or standard values were reached or exceeded.
- (26) Usability-Test, Usability Validation of AESCULAP® CeSPACE® 3D Cages, Tübingen, 2019. The usability of the AESCULAP® 3D Cage System CeSPACE® 3D was tested in April 2019, in a cadaver workshop with six independent test persons as intended users (surgeons specialized in spinal surgery or comparable fields). Parameters such as implant visibility under x-ray control, mechanical stability of the implant/instrument interface and implant surface evaluation in terms of tissue injury risk were tested among others. Acceptance criteria were fulfilled for all the above-mentioned parameters. All test users confirmed the absence of critical features that must be improved prior to clinical use. During the test, the x-ray visibility of the cages was particularly positively assessed.
- (27) Rehnitz, Christoph, PD Dr. med. Radiological image evaluation of AESCULAP® interbody fusion devices, Heidelberg, 2019. CT and X-ray visualization of different AESCULAP® interbody fusion cages (full titanium, porous Ti6Al4V and PLASMAPORE<sup>XP</sup> cages) was tested in a cadaver setup. A radiologist evaluated the implant visibility and the presence of artefacts that may limit the visualization of adjacent structures. Visualization and assessment of implant position was achieved in X-ray and CT for all tested cages. Minor artefacts were visible in CT reconstructions in the surrounding of porous Ti6Al4V and full titanium implants. Porous Ti6Al4V implants showed slightly fewer artefacts in CT in comparison to full titanium implants. The minor artefacts observed did not limit the assessment of the surrounding anatomical structures.

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[www.facebook.com/bbraunat](https://www.facebook.com/bbraunat) | [www.instagram.com/bbraun\\_at](https://www.instagram.com/bbraun_at) | [www.linkedin.com/company/bbraun-group](https://www.linkedin.com/company/bbraun-group)

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B. Braun Austria GmbH | Otto Braun-Straße 3-5 | 2344 Maria Enzersdorf
Hersteller nach MDR ist Aesculap AG, Am Aesculap-Platz, 78532 Tuttlingen, Deutschland